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Medimir – Clinical Decision Support System

Version: v2.0

UDI-DI: (01)4260756440109

UDI-PI:

- iOS: (8012) I020
- Android: (8012) A020

Date of Manufacturing: 2024.11.27

M E D I M I R

Contact:

support@medimir.ai

www.medimir.ai

Medimir Systems GmbH

Table of Contents - Inhaltsverzeichnis

Table of Contents - Inhaltsverzeichnis	2
English - Instructions for Use	4
Home Screen	4
Carousels	4
Navigation Bar	4
More Button	4
Info Menu	5
Tools Screen	5
List of Tools	5
Search	5
Patient List	6
App Settings	6
User account	7
User Menu	7
Automated Logout	7
Clinical Decision Support Tools	7
Notes, Bookmarks	8
Share, Print, Generate PDF	8
Basic Patient Data	9
Medication Plan	9
Laboratory Values	9
Management Board	10
Prediction Board	10
Recommendation Grades	10
Warnings and Precautions	11
Manufacturer, Label	12
Medical Device Information	13
Symbols	14

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English - Instructions for Use

Home Screen

The **Home screen** is the default view of the app and provides a quick overview and navigation to the various functions of the app.

The Home screen is organized by various **carousels**.

Access **specific areas** of the app by swiping through the carousels and selecting an area.

Tap an icon to access a resource directly.

View all available icons within a carousel by tapping "**Show all**".

Carousels

In the "**Newly released tools**" carousel, you find the tools recently released in the CDS system, each with a brief description.

The "**Top Tools Worldwide**" carousel gives you an overview of which tools are used most frequently by your colleagues.

The "**My Last Seen Tools**" carousel is a list of your personal top tools.

"**My Last Seen Patients**" lists the patient files you have worked on most recently.

You can tap on "**Show all**" at any time to see a complete list of the respective area.

Navigation Bar

The Navigation bar at the bottom of each screen provides quick and easy access to the main areas of the app via icons:

- (1) **Home**: Quick navigation back to home screen
- (2) **Patients**: List of all your patients
- (3) **Tools**: List of all interactive tools in the app and search function
- (4) **Chat**: Communicate with your colleagues via the chat function
- (5) **More**: Access Notes, Bookmarks and a History function

More Button

The More button of the bottom navigation bar gives you access to the following functions:

Search: Quickly navigate to the Search screen

Notes: Access your personal notes

Bookmarks: Manage and view your bookmarks

History: Easily navigate to previously viewed screens

Info Menu

Access additional functions through the Info menu in the top right-hand corner:

News: Read the latest news.

Instructions for Use: This introduction is a user manual that informs you about all the functions of the app.

Trainings: Watch the training videos on how to use the app correctly before using it.

App Settings: Choose your data backup, font size and update options and agree to data collection of your usage of the tools.

Technical support: Technical and content-related support requests or comments can be sent here directly to the developers.

Recommend: Share this resource via email, SMS or social media.

Rate the app: You can rate the app here.

Legal information: Information on data protection, disclaimer, imprint, copyright, etc.

Tools Screen

You can reach the **Tools screen** by tapping the Tools icon in the bottom navigation bar.

The **Tools screen** gives you access to a **Search** function (index, full text, filtered search) and the **List of Tools**, sorted by medical areas.

List of Tools

After choosing a medical area on the “**Tools Screen**”, tap on “All tools”, “CDS tools”, “Scores” or “Calculators” to narrow down the tools that will be displayed in an alphabetical order. You can access a specific tool by tapping on the respective title.

Search

Call up the **search function** by swiping down on the Home screen or tapping on the Tools icon in the lower navigation bar. You can then enter your search term into the search input field.

The **Search function** enables an “Index term”, “Full text” or “Filtered” search across the entire content of this app and your notes.

Enter your **search terms** in the search field and/or use the filter function in the top right-hand corner to refine your search.

- **Full text search:** performs a free full text search of the entire content of the app.
- **Index search:** searches only pre-indexed entries in the app.

- **Filtered search:** refines your search with specific filters.

Patient List

The **Patient list** can be accessed via the Patient button in the bottom navigation bar of the app.

All the patients you have created are listed in alphabetical order in this **Patient list**, either in a list view or in an icon view. Change the view of this patient list by tapping the List/icon view button.

Change the order of the patient list into listed „by date/time of admission“ by clicking on the button next to „Patients“.

At the top of the patient list there is a **search bar** to quickly search for a specific patient in the patient list.

You can **add a new patient** by clicking on the (“+”) plus-button in the right upper corner. This (“+”) plus-button opens a popup which requests for some minimum **basic patient data** (first name, last name) to create a new patient entry in the list. Any tool that is filled out and used for a particular patient, can now be saved into this patient file.

You can work on your patient files, use tools, save the results and manage all current patients within the „**Current**“ patients section.

By swiping the patient entry (in the patient list) to the left and tapping the „**Archive**“ button, the complete patient file can be moved into the „Archived“ section. The file can be retransferred into the „**Current**“ patient section any time.

A patient file can also be deleted by swiping the patient entry to the left and tapping the „**Delete**“ button.

App Settings

You can manage your **app settings** on the “**Settings**” page, accessible via the info menu.

Activate “**Data Backup**” so that your personal notes and bookmarks are not lost.

Select how you want to **update** the content. We recommend an **automatic update** via your mobile network.

Select the **language** of the app using the language switch.

Allow us to collect **analytics** data about your use of the tools to improve this app.

Change the **decimal separator** in the app tools according to your local standard.

Change the **font size** throughout the app.

User account

When **entering the CDS system** for the first time, the user is asked for their user account details. If available, the user can log in, if not available, the user will be asked to **register** and **set up a free Medimir user account**.

For reasons of **data security** and **data protection**, a **password** and a **four-digit PIN** must be created in addition to the user ID of the user account (your e-mail address). This data is required to be able to access the CDS system again later.

During the registration a range of **data about the user's professional activity** is collected. Please note that this software product is a **regulated medical device** and we must ensure that only healthcare professionals have access to this system. We therefore rely on the **accuracy of your selfdeclaration** that you are or will be a healthcare professional.

User Menu

Access your **personal user account** data within the app via the icon in the top right-hand corner of the app. Here you can **reset your password** or **PIN** or log out.

You can also **change or correct your personal data** that you entered when you first accessed the Medimir CDS system.

We use some of the user account data to **evaluate the usage of the system**, but only if you have authorized us to do so (see settings).

Automated Logout

If no interaction is registered for one minute while using the app, the **app will be locked**. A countdown timer and a warning are displayed 15 seconds before the app is locked so that the user can interrupt the lock. The app can be **unlocked** at any time using the **four-digit PIN**.

PIN and password can be reset or changed in the **personal user account area** within the app via the icon in the top righthand corner of the app.

Clinical Decision Support Tools

Clinical Decision Support tools (**CDS tools**) are software instruments that **guide you** step-by-step (and in the correct order) through **fulfilling guideline recommendations** in a specific patient case.

After **entering requested input data** about characteristics of the patient guideline-based **recommendations** for the management of the patient are **demonstrated to the user**.

Please ensure that **all required input** is entered to get correct results. The data items need to be selected from top to bottom and left to right.

On the left side of each recommendation, there is an arrow, that demonstrates the **level of recommendation**. Get additional information about data fields by tapping the **info buttons**.

Most questions of the CDS tools are set up in a simple **yes/no format**. Others require **selection of options** from a **dropdown menu** and some need manual input of **numbers** or text. Buttons can be selected by tapping on them, colors indicate the selection. Buttons can be deselected by tapping on them again.

When a **dropdown data item is clicked**, a list of options is shown on the screen. Tapping on one of the options will select it. In order to manually enter values on a free text entry data item (text or numerical value), tap the field and the device's system keyboard will appear.

Please ensure that **all required input** is entered into a tool to **receive correct output results**.

Depending on the answer to some of the questions, follow-up questions may be shown dynamically. As the questions are answered and decision rules are executed, subsequent sections of the tools will be filled with management recommendations dynamically. You can accept the recommendations by tapping "yes" or decline by tapping "no".

For **security reasons** the medical data on the device is **encrypted** and can be stored under a patient name in a cloud patient database. Only the user has access to these data, no one else. It is the user's responsibility to make sure the app and the device are not left unlocked and accessible for unauthorized persons to ensure confidentiality of the patient data.

In case the app user encounters any issues **with the system** components of the app, the user is obliged to **inform the manufacturer** in a timely manner.

Notes, Bookmarks

By clicking on the **three dots** in the title bar of a tool you can access additional functions:

- Create a **note** for the tool.
- Set a **bookmark** for the tool so that you can access it quickly.

By activating „**Data Backup**“ in the app settings your notes and bookmarks are automatically saved in the cloud of the CDS system.

You can access previously created notes and bookmarks via the **“More” icon in the bottom navigation bar**. You can also search the notes using the general search function.

Share, Print, Generate PDF

If you click on the **three dots** in the upper title bar of each tool, you can also:

- Read the **tool profile** or **tool information**.
- **Share** a link to a specific section of the tool, e.g. via AirDrop, Mail, LinkedIn, X.
- **Print** results or ...
- **Create a PDF** file of the respective section.

Basic Patient Data

Once a new patient has been added to the system (see patient list), the following **basic patient data** can be entered into the new patient file: First name, last name, gender, date of birth, height, weight, e-mail address and cell phone number.

Medication Plan

The **medication plan** contains all the medication prescribed for and taken by a patient. In the medication plan the medication is divided into continuous medication, ondemand, paused and discontinued medication.

Import medication plan

A patient's medication plan can be imported in full into this section of the app by scanning the QR code. Please tap on the QR code symbol at the top right and scan the code.

Add medication

A new medication can be added to the plan by tapping on the "Add line" drop-down menu and tapping on the "Add medication" option and then selecting it from an existing list based on ATC DDD codes. The user can then fill in all the relevant details in the medication form.

Discontinue, pause, resume, delete

Ongoing medication can be paused or discontinued as needed. Paused medication can be resumed or discontinued.

Laboratory Values

Laboratory values are entered manually into the laboratory values form.

The CDS system is already compatible (interoperable) with other information systems, so that the laboratory values in a hospital version of the system can also be exchanged automatically with other systems at a later date.

The laboratory values are divided into expandable and collapsible **category groups**.

Patient-specific laboratory values that are entered in this form are passed on to other tools of the CDS system. Other tools are then pre-filled with this data.

Management Board

The **management board** displays all the tools that have been used for a particular patient along with the recommendations of the tools after the user has finished entering inputs these tools.

Recommendations are labeled as **diagnostic** or **treatment** recommendations.

The dependencies between the recommendations are represented by nested indentations.

Prediction Board

The **prediction board** displays all calculated medical predictions resulting from prediction tools used for the respective patient. The displayed tools can be opened from this prediction board.

Recommendation Grades

Recommendation data items of the CDS tools are additionally equipped with a symbol that demonstrates the **grade (strength) of a recommendation**. Clicking on the symbol shows further information.

The following **three different levels of recommendation grades** are differentiated:

<green arrow> Strongly recommended, shall be considered: stands for „Evidence and/or general agreement that a given treatment or procedure is beneficial, useful, effective.

< green-bordered arrow> Recommended, should or may be considered: stands for „Conflicting evidence and/or a divergence of opinion about the usefulness/efficacy of the given treatment or procedure.

<black-bordered arrow> Not recommended: stands for Evidence or general agreement that a given treatment or procedure is NOT useful/effective, and in some cases may be harmful."

By clicking on the arrow symbol, the underlying narrative statement is shown and, if applicable, the exact grade of recommendation, the level of evidence and additional information about the source are also demonstrated.

Warnings and Precautions

	<p>For professional use only! Not a substitute for professional medical judgment! Verify information with primary sources!</p>
	<p>This Clinical Decision Support System is intended to support clinical decision making at the point of care. Although the CDS Tools are based on published guidelines and validated before release, it remains ultimately the user's responsibility to apply clinical reasoning when using the tools and to critically assess all recommendations before performing any action.</p>
	<p>The appropriate recommendation is shown only when all visible data items are activated. The user must answer all questions and activate the buttons accordingly to get the appropriate recommendation.</p>
	<p>The CDS tools can be used flexibly and do not force the user to perform actions in a specific order. However, for optimal use, it is recommended that the displayed columns are edited from left to right and each column from top to bottom.</p>
	<p>Incident Reporting: If you experience any serious incident that has occurred in relation to the device, please report it to us immediately via the Feedback form within this system and to the competent authority of the Member State in which you are established.</p>

Manufacturer, Label

Medimir - Clinical Decision Support System

Version: v2.0

Date of manufacture (release date): 2024.11.27

UDI-DI: (01)4260756440109

UDI-PI: (8012) IO20

Manufacturer:

Börm Bruckmeier Verlag GmbH,
Emil-Geis-Str. 4, 82031 Gruenwald, Germany

Product Owner, Copyright:

© 2024, Medimir Systems GmbH,
Emil-Geis-Str. 4, 82031 Gruenwald, Germany

Contact:

support@medimir.ai

www.medimir.ai

These **instructions for use** are only provided in electronic form in accordance with Article 6(1) of Regulation (EU) No 207/2012.

eIFU Access: Current and historical versions of the **instructions for use** are available at www.medimir.ai.

For **support**, please contact support@media4u.com.

Medical Device Information

Name: Medimir - Clinical Decision Support System

Description: Mobile application providing clinical decision support for healthcare professionals

Medical Device Status: This App is a medical device as defined by the Medical Device Regulation (MDR) 2017/745, Article 2(1). Its medical purpose is to support diagnosis, prevention, monitoring, prediction, prognosis, treatment or alleviation of disease.

Intended Purpose: As a medical device, this App is designed to support healthcare professionals in managing clinical knowledge and making guideline-compliant clinical decisions at the point of care, taking into account up-to-date patient data.

Intended User: This medical device software is intended to be used by healthcare professionals.

Intended Location of Use: This medical device app is intended for use in hospitals and doctors' practices.

Medical Indications and Intended Patient Population: As a multidisciplinary medical device, there are no restrictions on specific indications for the use of this app in relation to anatomical structures, diseases, injuries, or specific medical fields. There is no restriction regarding the intended patient population.

Contraindications: None

Clinical Benefits: quick access to current medical information

Residual risks: potential for outdated information if not updated regularly

Classification according to Regulation (EU) 2017/745 on medical devices: This software/app system is classified as a Class IIa medical device under Rule 11 of Annex VIII.

Software Safety Classification according to IEC 62304: This medical device software system is classified as Class B software.

Software Requirements: iOS 12.0 or later / Android 8.0 or later, internet connection required for updates

Installation: download from App Store or Google Play Store

Updates: automatic or manual via app stores

Security: regular updates for cybersecurity, no patient data stored or transmitted

Data and privacy: no personal data collected or stored

Symbols

	<p>This symbol indicates a warning or important note, which is written next to it.</p>
	<p>This symbol indicates the date of manufacture of the device (medical device).</p>
	<p>This symbol indicates the manufacturer of the device (medical device).</p>
	<p>This symbol indicates the manufacturer's declaration that the product meets EU standards for health, safety, and environmental protection. The four digit identification number indicates the Notified Body involved in the conformity assessment procedure.</p>

Abbreviations

Abbreviation	Meaning
BB	Börm Bruckmeier
CDS	Clinical Decision Support
CDSS	Clinical Decision Support System
CDST	Clinical Decision Support Tool
POC	Point of Care

Indication-specific clinical abbreviations are used in many of the CDS tools of this CDS system. These abbreviations are always provided within the tool in a separate, easily accessible list of abbreviations.

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